

# Transducers

C5-2Q



2 - 5 MHz  
Curved Array

C5-1Q



1 - 5 MHz  
Single Crystal  
Curved Array

MC8-4Q




4 - 8 MHz  
Micro - Convex Array

MC9-3TQ



3 - 9 MHz  
Micro - convex Array

C6-2MQ



2 - 6 MHz  
3D/4D Mechanical  
Curved Array

L12-5Q



5 - 12 MHz  
Linear Array

L17-7HQ



7 - 17 MHz  
High Frequency  
Linear Array

L17-7SQ



7 - 17 MHz  
High Frequency Compact  
Linear Array

E10-3BQ



3 - 10 MHz  
Endocavity  
Curved Array ( 200° )

E10-3HQ



3 - 10 MHz  
Endocavity  
Curved Array ( 200° )

P5-1Q



1 - 5 MHz  
Phased Array


P7-3Q



3 - 7 MHz  
Pediatric Phased Array

Being Productive  
Being Trust

Acclarix™ LX3  
Diagnostic Ultrasound System



Specializes in the research and development, manufacturing, sales and service of technologically advanced medical devices, which products cover medical imaging, cardiovascular and peripheral minimally invasive surgery, anesthesia, intensive care and other departments, and vaccine syringes (such as the luer-lock type for Covid-19).

911 SW 3RD ST APT1, MIAMI, FL 33130-2227

+786-9201598

sales@globalmedica.us

EDAN  
A world of potential






# Stunning Clarity

Compacted with innovative Edan TAI technology and multiple imaging processing technologies, Acclarix LX3 could perfectly display ultrasound image in different modes, assisting sonographer to make more precise diagnosis.

- High fidelity, high-channel count architecture results in superb detail resolution, particularly at depth
- Tissue Adaptive Imaging (TAI) continuously and automatically optimizes imaging allowing more focus on the patient
  - In B-mode, TAI fine tunes multiple parameters to provide the best possible image quality
  - In Doppler, TAI automatically adjusts for flow state, providing improved continuity, border detection and filling of color



Fetal Face




Fetal Heart




MSK



Color Doppler Spleen



TVD



Anatomic M

# New Generation of Capability and Flexibility

Integrated with compact appearance, cutting-edge technology and intelligent workflow, Acclarix LX3 is to provide more valuable innovation. Its excellent image performance plays an important role in the prevention and early detection in primary medical care, enriching the diagnostic application and enhancing the clinical value.

# The Virtue of Value

The Acclarix LX3 delivers unmatched value and performance across a broad range of applications:

- Abdomen
- Cardiac Screening
- Small Parts
- OB/GYN
- Vascular
- MSK
- Difficult-to-Image
- Breast
- Nerve

# Distinctive Design



Considerate gel warmer design for better experience



Customizable touch screen user panel to simplify diagnosis



5 active transducer ports to maximize clinical application



Keyboard with dedicated swivel for easy folding