



Pressure Connecting Tubes



Specializes in the research and development, manufacturing, sales and service of technologically advanced medical devices, which products cover medical imaging, cardiovascular and peripheral minimally invasive surgery, anesthesia, intensive care and other departments, and vaccine syringes (such as the luer-lock type for Covid-19).

911 SW 3RD ST APT1, MIAMI, FL 33130-2227

+786-9201598

sales@globalmedica.us

antmed
Partner For Life

Pressure Connecting Tubes

- Ensure Patient Safety - DEHP Free, comprehensive clinical validation and stringent bio-compatibility testing
- Enhance Operational Efficiency - Excellent in temperature and pressure resistant, the deformation of tube is smaller, yields great effect for angiography
- A variety of styles of tubing to best fit your needs - Lengthened, coiled, straight, T-tube, one-way checkvalve and backflow prevention device. Customization requests can be accommodated



• Tubes for CT (straight & coiled)



• Tubes for MR (Y & T type)



• Tubes for CathLab (600 psi)



• Tubes for CathLab (1200 psi)

Connecting Tube Type	Length	Pressure Range MPa(psi)	Specifications
CT Connecting Tube	5cm ~ 250cm	2.4(350)	Y & T straight, coiled tube with checkvalve(s), stopcock & spike
MR Connecting Tube	5cm ~ 250cm	2.4(350)	Y & T straight, coiled tube with checkvalve(s), stopcock & spike
CathLab Connecting Tube	5cm ~ 250cm	4.2(600)	Straight line with rotating male Luer-Lock
	5cm ~ 250cm	8.3(1200)	Braided line with rotating male Luer-Lock



Specializes in the research and development, manufacturing, sales and service of technologically advanced medical devices, which products cover medical imaging, cardiovascular and peripheral minimally invasive surgery, anesthesia, intensive care and other departments, and vaccine syringes (such as the luer-lock type for Covid-19).

911 SW 3RD ST APT1, MIAMI, FL 33130-2227

+786-9201598

sales@globalmedica.us