






High Pressure IV Catheter

- 350PSI
- Designed for CT and MR power injection of contrast

Specifications and Colors

Specification	Color	Maximum pressure
24G×0.7mm×19mm(24G×3/4")	Yellow 	350PSI
22G×0.9mm×25mm(22G×1")	Deep blue 	350PSI
20G×1.1mm×29mm(20G×1 1/8")	Pink 	350PSI
18G×1.3mm×29mm(18G×1 1/8")	Deep green 	350PSI

Specializes in the research and development, manufacturing, sales and service of technologically advanced medical devices, which products cover medical imaging, cardiovascular and peripheral minimally invasive surgery, anesthesia, intensive care and other departments, and vaccine syringes (such as the luer-lock type for Covid-19).

911 SW 3RD ST APTL MIAMI, FL 33130-2227

+786-9201598

sales@globalmedica.us